



DATASHEET
EOC-SI-2400
High Performance Ultra Thin Fiber Spectrometer

Features:

- The minim, ultra-thin, ultra-light weight of a minim;
- Fiber input and signal output are on the same side, which is very easy to integrate.
- Detector: linear CMOS
- Detector pixel: 2048 pixels
- Ultra-low noise CCD signal processing circuit
- Maximum spectral range: 200-1180 nm (depending on requirements)
- Spectral resolution: 0.1-2 nm (depending on spectral range, slit width)
- Optical Design: cross-C-T
- Built-in pulsed xenon lamp driver
- Casual time: 0.1ms-130s
- Vinner power supply: DC 5V ±10% or USB power supply
- A16 bit, 2MHz ADC
- Incident Interface: SMA905 or free space
- Supported data output interface: USB 2.0 (High speed) or UART
- 20-pin dual-row programmable external expansion interface

Description:

After 20 years in the development of optical fiber spectrometers, a new generation of high-performance ultra-thin fiber spectrometer, the EOC-SI-2400, is now available. It has built-in pulsed xenon lamp drive circuits, adopts high sensitivity linear CMOS, and specially customized ultra-low noise CMOS signal processing circuits. The noise of the sensor is lower, the signal-to-noise ratio is obtained (about twice higher than that of similar competitors), and the measurement reliability of the EOC-SI2400 is improved. The measurement results do not change with the ambient temperature, which is the best level in the industry.

EOC-SI-2400 can receive SMA905 fiber input light or free space light, and output spectral data measured through USB2.0 or UART ports.

EOC-SI-2400 only requires a 5V DC power supply or USB power supply, which is very easy to integrate.

| Model | Detector pixel | Sensor Cooled |
|-------------|----------------|---------------|
| EOC-SI-2400 | 2048pixel | No |



Application:

- Fast and ultra-micro spectrophotometer; ;
- Environmental protection equipment (flue gas, water quality);
- Reflection and transmission spectrum detection;
- Spectral analysis, radiation spectroscopic analysis, spectrophotometric analysis;
- Laser wavelength measurement.

1. Performance parameters:

| Detector | |
|------------------------------|---|
| Type | Linear CMOS |
| Detectable range | 200-1100 nm |
| Effective pixel | 2048×1 |
| Pixel dimension | 14μm×200μm |
| Full scale range | ~200 ke- |
| Sensitivity | 1300 V/(lx•s) |
| Dark noise | 0.4 mV rms |
| Optical Parameter | |
| Wavelength range | 200-1100 nm (Specific according to needs) |
| Optical resolution | 0.1-2 nm (Depends on slit, spectral range) |
| SNR | >450 : 1 |
| Dynamic range | 3000: 1 |
| Optical Configuration | |
| Optical Design | f/4 cross asymmetric C-T optical path |
| Incidence slit | 5、 10、 25、 50、 100、 150、 200 μm Optional, other sizes can be customized |
| Incident Interface | SMA905 Optical fiber interface, free space |
| Electrical Parameter | |
| Integration time | 0.1 ms - 256 second |
| Interfaces | USB 2.0 |
| A/D conversion resolution | 16 bit |
| Supply voltage | DC4.5 to 5.5 V (type @5V) |
| Operating current | 250mA@Typ. |
| Operating current | -30°C to +70°C |
| Operating temperature | -25-50°C |
| Working humidity | < 90%RH (No condensation) |
| Physics Parameter | |
| Dimension | 85 × 61 × 24 mm |
| weight | 0.15 kg |

2 Mechanical Diagrams



3 Electrical Pin-out

Table 1 Electrical Characteristics

| Parameter | Min | Typ | Max | Unit |
|---|------|-----|-----|------|
| Power Supply | | | | |
| Operating voltage range | 4.5 | 5 | 5.5 | V |
| Operating current | | 250 | | mA |
| Logic Inputs(3.3V LVTTL, Five-volt tolerant) | | | | |
| High level input voltage | 1.7 | | 3.6 | V |
| Low level input voltage | -0.3 | | 1.0 | V |
| Logic Output(3.3V LVTTL) | | | | |
| High level output voltage | 2.4 | | 0.4 | V |
| Low level output voltage | | | | V |

The module is equipped with a 20-pin male angled box header(2x10, 2.00 mm pitch) and USB2.0 C type interface. The 20-pin connector is a Samtec part # STMM-110-02-L-D-RA connector. The mate to this is a Samtec part # TCSD-10-D-XX.XX-01-N.




Table 2 Electrical Pin-Out

| Pin# | Description | I/O | Function Description |
|------|----------------|--------|---|
| 1 | VCC | / | Power Supply, 5V±0.5, |
| 2 | GND | / | Ground |
| 3 | TX | Output | LVTTL Transmit signal |
| 4 | RX | Input | LVTTL Receive signal |
| 5 | Lamp Enable | Output | LVTTL output the lamp enable signal. |
| 6 | N.C. | / | / |
| 7 | Ext_trigger_in | Input | LVTTL input the trigger signal. Falling edge trigger collection |
| 8 | N.C. | / | / |
| 9 | SPI_SCK | Output | The SPI Clock signal for communications to other SPI peripherals |
| 10 | SPI_MOSI | Output | The SPI Master Out Slave In (MOSI) signal for communications to other SPI peripherals |

| | | | |
|----|----------|---------------|---|
| 11 | SPI_MISO | Input | The SPI Master In Slave Out (MISO) signal for communications to other SPI peripherals |
| 12 | SPI_CS | Output | The SPI Chip/Device Select signal for communications to other SPI peripherals |
| 13 | GPIO0 | Input /Output | General Purpose Software Programmable Digital Inputs/Outputs, LVTTTL Logic. |
| 14 | GPIO1 | Input /Output | General Purpose Software Programmable Digital Inputs/Outputs, LVTTTL Logic. |
| 15 | GPIO2 | Input /Output | General Purpose Software Programmable Digital Inputs/Outputs, LVTTTL Logic. |
| 16 | GPIO3 | Input /Output | General Purpose Software Programmable Digital Inputs/Outputs, LVTTTL Logic. |
| 17 | GPIO4 | Input /Output | General Purpose Software Programmable Digital Inputs/Outputs, LVTTTL Logic. |
| 18 | GPIO5 | Input /Output | General Purpose Software Programmable Digital Inputs/Outputs, LVTTTL Logic. |
| 19 | GPIO6 | Input /Output | General Purpose Software Programmable Digital Inputs/Outputs, LVTTTL Logic. |
| 20 | GPIO7 | Input /Output | General Purpose Software Programmable Digital Inputs/Outputs, LVTTTL Logic. |

4 Spectrogram example

