



**DATASHEET  
EOC-SI-R8300P  
Raman Microscope (Pro)**

## Description

EOC-SI-R8300P is a miniature Raman microscope integrating benefits of microscope and Raman spectrometer into one instrument. It becomes possible to see micro areas of samples on the computer screen with just a single mouse click. When the sample is visualized in an accurate position, the observer can scan Raman spectrum under various surface conditions, and synchronous Mapping can be displayed intuitively on the screen in a click. As a result, it takes great convenience to detect micro areas of samples. Combine unique patented conjugate focusing (true confocal) system with accurate image processing algorithm, and it enables a very small sample areas to be analyzed, as well as it requires minimal operator training and maintenance, yet resulting in uniform result not just spectrum.

EOC-SI-R8300P is equipped with tailor-made objective, and laser spot on the sample becomes very close to diffraction limit, then focal information can be displayed in accurate and intuitive on the screen with 3-megapixel/5-megapixel camera. This configuration improves Raman spectral quality for overcoming the limitations of Raman systems where the focal plane for Raman signal collection is slightly above or below the imaging plane.

EOC-SI-R8300P works very stable with no moving components of optical path switch, hence it avoids loss off optical path while imaging being formed, and it gains optimized signal for separating imaging formed from Raman signal collection.

## Feature:

- Full-automated, auto-focusing, auto-scan
- Ultra-high resolution  $1\text{cm}^{-1}$ .
- Ultra-high sensitivity  $>6000:1$
- True confocal, accurate Raman mapping
- Unique software controlled to switch optical path
- Ultra-high stability
- Excellent performance
- Fast positioning, quick locate focal position
- High quality objective, micro spot
- 3-mega/5-mega pixel camera, crisp clear images
- Excitation wavelength(Optional): 532,633, 785,830,1064
- High-performance spectrometer configured
- USB2.0 in direct connect with PC

## Application:

- Nano particles and new materials
- Science research Institutions
- Bioscience
- Forensic identification
- Material science
- Medical immunology analysis
- Agriculture and food accreditation
- Gemstones & minerals identification



Fig. 1 EOC-SI-R8300Pro Structure indication diagram

Table 1 EOC-SI-R8300 Pro Product Selection

| Models                     | Focus Length | Excitation Wavelength/nm | Excitation Power/mW | Max. Wavenumber Range | Min. Resolution/cm <sup>-1</sup> |
|----------------------------|--------------|--------------------------|---------------------|-----------------------|----------------------------------|
| EOC-SI-R8300Pro<br>-FL 210 | 210mm        | 532                      | 100                 | 200-3500              | 2.2                              |
|                            |              | 633/638                  | 80                  | 200-3300              | 2.2                              |
|                            |              | 785                      | 350                 | 200-3500              | 2.5                              |
|                            |              | 1064                     | 500                 | 200-2500              | 6.2                              |
| EOC-SI-R8300Pro<br>-FL 350 | 350mm        | 532                      | 100                 | 200-3700              | 1.4                              |
|                            |              | 633/638                  | 80                  | 200-3500              | 1.4                              |
|                            |              | 785                      | 350                 | 200-3500              | 2.1                              |
|                            |              | 1064                     | 500                 | 200-2500              | 5.1                              |
| EOC-SI-R8300Pro<br>-FL 510 | 510mm        | 532                      | 100                 | 200-3700              | 0.9                              |
|                            |              | 633/638                  | 80                  | 200-3500              | 0.9                              |
|                            |              | 785                      | 350                 | 200-3500              | 1.4                              |
|                            |              | 1064                     | 500                 | 200-2500              | 3.6                              |
| EOC-SI-R8300Pro<br>-FL 760 | 760mm        | 532                      | 100                 | 200-3700              | 0.5                              |
|                            |              | 633/638                  | 80                  | 200-3500              | 0.5                              |
|                            |              | 785                      | 350                 | 200-3500              | 1.0                              |

|   |  |      |     |          |     |
|---|--|------|-----|----------|-----|
|   |  | 1064 | 500 | 200-2500 | 2.7 |
| EOC-SI-R8300Pro-LT: Cooled down to -30°C, ultra-long integration time (Max. Time can reach 1.3h).<br>EOC-SI-R8300Pro-SCM: Te-cooled SCMOS detector.<br>EOC-SI-R8300Pro-BS: Basic type.<br>EOC-SI-R8300Pro-AF: Auto-focus.<br>EOC-SI-R8300Pro-MP: Mapping, and auto-focus. |  |      |     |          |     |

Naming example:

EOC-SI-R8300Pro-AF-LT-FL350-532+633: auto focus, long integration time, focus length of 350mm, dual excitation wavelength: 532nm and 633nm respectively

EOC-SI-R8300Pro-MP-SCM-FL760-532+633+1064: scanning imaging, SCMOS detector, focus length is 760mm, excitation wavelength is three wavelengths: 532nm, 633nm and 1064nm respectively

## 1. Specification

| EOC-SI-R8300Pro                            |  |
|--|--|
| Excitation wavelength                      | Refer to table 1.  |
| Spectral resolution                        | Refer to table 1.  |
| Spectral range                             | 250-2700, 200-3500, 200-4300 cm-1<br>(available in customer wavelengths range down to 50 cm-1)           |
| Maximum laser output                       | 500mW (Max. 100mW for 532nm)   |
| Spectral Stability                         | $\sigma/\mu < 0.5\%$ (COT 8 hours)   |
| Thermal stability                          | Spectral shift $\leq 1$ cm-1 (10-40 °C)  |
| SN ratio                                   | >6000:1  |
| Detector                                   | TE cooled, semiconductor laser, 2048*64 pixel, back-thinned, IR enhanced CCD<br>InGaAS cooled for 1064nm |
| wavelength range detected                  | 200nm-1100nm   |
| Pixel size                                 | 14 $\mu\text{m}$ * 14 $\mu\text{m}$  |
| Dynamic range                              | 13000:1  |
| Laser center wavelength                    | 785nm (+/-0.5nm)   |
| Microscope camera                          | 3-megapixel /5-megapixel camera  |
| focusing                                   | True confocal  |
| Laser output                               | >550mW (software adjustable)   |
| laser spot diameter                        | >1 $\mu\text{m}$   |
| Laser stability                            | $\sigma/\mu < \pm 0.2\%$   |
| Laser linewidth                            | 0.08 nm  |
| Connectivity                               | USB2.0   |
| Electrical controlled X,Y axis 2D platform |  |

|                             |                     |
|-----------------------------|---------------------|
| moving range                | 5 X 5 cm            |
| moving resolution           | 0.1 $\mu$ m         |
| positioning accuracy        | 1 $\mu$ m           |
| Scan speed                  | 20mm/s              |
| Z axis (automated focusing) |                     |
| focusing accuracy           | $\leq \pm 0.2\mu$ m |
| Max. range                  | 20mm                |
| focusing speed              | Less than 10 s      |

## 2. Optical performance

### 2.1 Spectrum








Fig. 1 EOC-SI-R8300Pro collect acetonitrile spectra






**Fig.2 EOC-SI-R8300Pro experiment: Left picture is sample, and right picture is Sers Raman spectra**



**Fig 3 EOC-SI-R8300Pro Measure Si Raman spectra (500mW, integration time: 1S)**



**Fig 4 EOC-SI-R8300Pro measure alcohol spectra (500mW, integration time:1S)**




Fig 5 EOC-SI-R8300Pro operation interface




Fig 6 EOC-SI-R8300Pro operation interface

## 2.2 Raman resolution

### 2.2.1 Tylenol Raman spectra




Fig 7 Tylenol spectra shows clear 1610/1615  $\text{cm}^{-1}$  vibration peak

### 2.2.2 Petrol Raman spectra

Power output: 200 mW  
Integration time: 10 s  
Filter level 1

Raman spectra of 93# petrol showed the resolution condition in the long wavelength region.




Fig 8 93# petrol Raman spectra, 723/732/742 $\text{cm}^{-1}$  spectral shift is clearly recognized

### 3. Reliability

Fig 3.1, Fig 3.2 temperature stability is measured by EOC-SI-R8300, kept stable above an hour for each temperature node ranging between 5-40°C. Sample measured is acetonitrile, wavenumbers shift  $\leq 1\text{cm}^{-1}$ (Fig 3.1), peak top intensity change  $< 10\%$  (Fig 3.2)






Fig. 3.1 Wavenumber shift results testing from 5 °C to 40 °C of five ATR2000 portable Raman spectrometers




Fig. 3.2 Intensity variation testing from 5 °C to 40 °C of five ATR2000 portable Raman spectrometers