

# Electro Optical Components, Inc.

5460 Skylane Boulevard, Santa Rosa, CA 95403 Toll Free: 855-EOC-6300 www.eoc-inc.com | info@eoc-inc.com



## DATASHEET EOC-SI-R3000 High Sensitivity & High Resolution Portable Raman Spectrometer

#### Feature:

- Ultra-high sensitivity FFT-CCD TE-cooled;
- low noise circuit;
- Powerful embedded software;
- Fluorescent background eliminate;
- Peak finding and display;
- ➤ 10.1-Inch LCD;
- Win 10 operation system;
- > 11.6-inch capacitive touch screen, multitouch;
- USB 2.0;
- User friendly human-machine interface;
- Battery life> 3h;
- Remote control via LAN;
- IP67 case;

#### Application:

- Biological science
- Pharmaceutical engineering
- Forensic analysis
- Agriculture and food safety
- Gemstone
- Environmental science

### **Description:**

EOC-SI-R3000 portable Raman spectrometer is suitable for field operation. The outstanding reliability makes the detection result much accurate.

The excellent low stray light conditions that enable the spectrometer has a wide range of application, especially in public safety, food safety, and pharmaceutical engineering.

The multi-function software facilitated the spectral analysis process in application. The remote experiment through internet access makes the test item much easier.

ltem	Wavelength	resolution
	range (cm <sup>-1</sup> )	(cm <sup>-1</sup> ) *
EOC-SI-R3000-27	150-2600	5
EOC-SI-R3000-35	200-3500	5
EOC-SI-R3000-43	200-4300	6

Remark:

- Measuring method is based on ASTM E2529-06;
- Available in custom design, resolution can be increased by around 1/3, resulting in lower sensitivity;



# Specifications

EOC-SI-R3000			
Interface	USB 2.0 and WIFI		
Operating system	Android		
Screen	11.6-inch capacitive touch screen, Multi-touch		
Battery life	>4 h		
Integration time	4ms - 120s		
Power voltage	DC 19V(+/-5%)		
Operating Temp	-10~40 °C		
Operating humidity	< 95%		
Dimension(L*W*H)	40×30×18 cm3		
Weight	7.5 Кg		
Reliability			
Spectral stability	σ / μ < 0.5% (COT 8 hours)		
Temp stability	Spectral shift $\leq$ 1 cm-1 (10-40 $$ $^{\circ}\mathrm{C}$ )		
Variation of intensity (in 5 ~ 40 $^{\circ}$ C)	<±5%		
Optical parameters			
Spectral range (cm-1)	150-2600	200-3500	200-4300
resolution (cm-1)	5	5	6
SNR	>3000:1 (918 cm-1 of Acetonitrile , 10s accumulation, 200mW)		
Entrance slit	50 μm		
Optical system	f/4 C-T crossed optical path		
focusing	98 mm for incidence and output		
Detector			
ltem	Ultra-high sensitivity, quick cooling CCD		
Detector cooled down to	-10 °C		
Detecting range	200-1100 nm		
Effective pixels	2048*64		
Dynamic range	50000: 1		
Pixel size	$14 \mum  imes 14 \mum$		
Full well capacity	300 Ке-		
Sensitivity	QE>40%, 6.5 μ V/e-		
Exciting Laser			

Central wavelength	785nm (+/-1nm)
FWHM	0.08 nm
Power output	≥500 mW
Power stability	σ/μ <±0.2%
Raman probe	
Operating distance	6 mm
Rayleigh scattering resistance	OD>8
Numerical Aperture	0.3
Aperture	7mm



Figure 1 EOC-SI-R3000