

Electro Optical Components, Inc.

5464 Skylane Boulevard, Suite D, Santa Rosa, CA 95403 Toll Free: 855-EOC-6300

www.eoc-inc.com | info@eoc-inc.com







MarginGuide for Tumor Diagnostics

Raman Spectroscopy Diagnostics ex-vivo to detect tumor margins in 3D using the thinnest fiber probe

Raman Needle Probe

penetrates in resected tissue ex-vivo to define oral cancer (SCC) margins for 1cm depth within 1s – measuring each 100ms the ratio of water/lipids bands in HW-Raman spectrum (as its growth is SCC specific biomarker)







